

Agpat3 Rabbit Polyclonal Antibody For Wb And Elisa - Q9NrZ7

Item Number: CM0013541

Introduction

Rabbit polyclonal antibody against human AGPAT3, validated in WB and ELISA with cross-reactivity to mouse. Targets 1-acylglycerol-3-phosphate acyltransferase gamma involved in phospholipid biosynthesis. Suitable for research on lipid metabolism.

[Learn More](#)

Parameter	Value
Product Name	AGPAT3 Rabbit pAb
Remarks/Alias	LPAAT3; LPLAT3; 1-AGPAT 3; LPAAT-GAMMA1; AGPAT3
Species	Human
Gene ID (Human)	56894
Gene ID	56894
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9NRZ7
Protein Weight	43kDa
Shipping	Ice bag