

Agr2 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O95994

Item Number: CM0030844



Introduction

AGR2 rabbit polyclonal antibody validated for Western blot, immunohistochemistry, immunofluorescence, and ELISA. Recognizes human, mouse, rat anterior gradient protein 2, a proto-oncogene involved in mucus synthesis and cell adhesion.

[Learn More](#)

Parameter	Value
Product Name	AGR2 Rabbit pAb
Remarks/Alias	AG2; AG-2; HPC8; GOB-4; HAG-2; RIFTD; XAG-2; PDIA17; HEL-S-116; AGR2
Species	Human
Gene ID (Human)	10551
Gene ID	10551
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-175 of human AGR2 (NP_006399.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95994
Protein Weight	20kDa
Shipping	Ice bag