

Ahnak Rabbit Pab For Ihc-P, If/Icc, Elisa - Q09666

Item Number: CM0024043



Introduction

High-quality rabbit polyclonal antibody against human AHNAK (Neuroblast differentiation-associated protein AHNAK), validated for IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for neuroscience and cell differentiation research.

[Learn More](#)

Parameter	Value
Product Name	AHNAK Rabbit pAb
Remarks/Alias	PM227; AHNAK1; AHNAKRS; AHNAK
Species	Human
Gene ID (Human)	79026
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-149 of human AHNAK (NP_001611.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q09666
Protein Weight	629kDa
Shipping	Ice bag