

Aif1/Iba1 Rabbit Monoclonal Antibody - 070200

Item Number: CM0002587

Introduction

Rabbit monoclonal antibody against mouse AIF1/IBA1. Validated for WB, IF, IHC-P, and ELISA, with cross-reactivity to human and rat. Suitable for studying allograft inflammatory factor 1 in immune cells and microglia.

[Learn More](#)

Parameter	Value
Product Name	AIF1/IBA1 Rabbit mAb
Remarks/Alias	G1; Iba1; AIF-1; D17H6S50E; AIF1/IBA1
Species	Mouse
Gene ID (Human)	199
Gene ID (Mouse)	11629
Immunogen	Synthetic Peptide, A synthetic peptide corresponding to a sequence within amino acids 1-100 of mouse AIF1/IBA1 (070200).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-F, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	070200(mouse), P55008(Human)
Protein Weight	17kDa
Shipping	Ice bag