

Aipl1 Rabbit Polyclonal Antibody (Pab) For Wb, If/Icc, Elisa - Q9Nzn9

Item Number: CM0013993



Introduction

High-quality AIPL1 rabbit polyclonal antibody for Western blot, immunofluorescence, and ELISA applications. Detects endogenous AIPL1 in human, mouse, and rat samples. Ideal for studying protein trafficking and Leber congenital amaurosis.

[Learn More](#)

Parameter	Value
Product Name	AIPL1 Rabbit pAb
Remarks/Alias	LCA4; AIPL2; AIPL1
Species	Human
Gene ID	23746
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-384 of human AIPL1 (NP_055151.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NZN9
Protein Weight	44kDa
Shipping	Ice bag