

Ak5 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Y6K8

Item Number: CM0024669



Introduction

High-quality AK5 rabbit polyclonal antibody validated for WB, IHC-P, ELISA with cross-reactivity to human, mouse, rat. Ideal for studying adenylate kinase 5 in brain tissue.

[Learn More](#)

Parameter	Value
Product Name	AK5 Rabbit pAb
Remarks/Alias	AK6; AK5
Species	Human
GeneID (Human)	26289
GeneID	26289
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 150-370 of human AK5 (NP_777283.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y6K8
Protein Weight	63kDa
Shipping	Ice bag