

Akt1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Ip - P31749

Item Number: CM0003939



Introduction

High-affinity AKT1 rabbit monoclonal antibody validated for WB, IHC-P, IP. Recognizes human, mouse, rat. Targets RAC-alpha serine/threonine kinase, a central regulator of cell survival, metabolism, and proliferation. Ideal for cancer and signaling research.

[Learn More](#)

Parameter	Value
Product Name	AKT1 Rabbit mAb
Remarks/Alias	AKT; PKB; RAC; PRKBA; PKB-ALPHA; RAC-ALPHA; AKT1
Species	Human
Gene ID (Human)	207
Gene ID	207
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 381-480 of human AKT1 (P31749).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P31749
Protein Weight	56kDa
Shipping	Ice bag