

Akt1 Rabbit Monoclonal Antibody (Ko Validated) For Wb, Ihc-P, If/Icc, Elisa - P31749

Item Number: CM0004818



Introduction

KO Validated AKT1 Rabbit Monoclonal Antibody specifically detects AKT1 in human, mouse, and rat samples. Ideal for Western Blot, IHC-P, IF/ICC, and ELISA applications in AKT signaling research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] AKT1 Rabbit mAb
Remarks/Alias	AKT; PKB; RAC; PRKBA; PKB-ALPHA; RAC-ALPHA; T1
Species	Human
Gene ID (Human)	207
Gene ID	207
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 380-480 of human AKT1 (NP_005154.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P31749
Protein Weight	56kDa
Shipping	Ice bag