

## Aldh3A2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P51648

Item Number: CM0022707



### Introduction

ALDH3A2 Rabbit Polyclonal Antibody validated for WB, IF/ICC, ELISA in Human, Mouse, Rat. Targets fatty aldehyde dehydrogenase, essential for lipid metabolism and sphingosine-1-phosphate degradation. Ideal for metabolic disease and Sjögren-Larsson syndrome research.

[Learn More](#)

| Parameter        | Value                                                                                                                   |
|------------------|-------------------------------------------------------------------------------------------------------------------------|
| Product Name     | ALDH3A2 Rabbit pAb                                                                                                      |
| Remarks/Alias    | SLS; FALDH; ALDH10; ALDH3A2                                                                                             |
| Species          | Human                                                                                                                   |
| GeneID (Human)   | 224                                                                                                                     |
| GeneID           | 224                                                                                                                     |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 250-320 of human ALDH3A2 (NP_001026976.1) |
| Source           | Rabbit                                                                                                                  |
| Category         | Polyclonal Antibodies                                                                                                   |
| Application      | WB, IF/ICC, ELISA                                                                                                       |
| Cross Reactivity | Human, Mouse, Rat                                                                                                       |
| SWISS            | P51648                                                                                                                  |
| Protein Weight   | 55kDa                                                                                                                   |
| Shipping         | Ice bag                                                                                                                 |