


Aldh9A1 Rabbit Pab - P49189

Item Number: CM0029877



Introduction

The ALDH9A1 Rabbit Polyclonal Antibody is validated for WB, IF/ICC, and ELISA, detecting human, mouse, and rat ALDH9A1. Ideal for aldehyde dehydrogenase research, this antibody targets an enzyme involved in biogenic amine and polyamine metabolism.

[Learn More](#)

Parameter	Value
Product Name	ALDH9A1 Rabbit pAb
Remarks/Alias	E3; ALDH4; ALDH7; ALDH9; TMABADH; TMABA-DH; TMABALDH; ALDH9A1
Species	Human
GeneID (Human)	223
GeneID	223
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human ALDH9A1 (NP_000687.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49189
Protein Weight	54kDa
Shipping	Ice bag