

Alk Rabbit Mab For Wb, Ihc-P, Elisa - Q9Um73

Item Number: CM0012088



Introduction

Recombinant rabbit monoclonal antibody against human ALK (CD246) for WB, IHC-P, ELISA. Detects anaplastic lymphoma kinase, a neuronal receptor tyrosine kinase crucial in nervous system development and cancer biology.

[Learn More](#)

Specification	Details
Product Name	ALK Rabbit mAb
Remarks/Alias	ALK1; CD246; NBLST3
Species	Human
Gene ID (Human)	238
Gene ID	238
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1093-1411 of human ALK (NP_004295.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	Q9UM73
Protein Weight	176 kDa
Shipping	Ice bag