


Alpha-1 Antichymotrypsin (Serpina3) Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P01011

Item Number: CM0028387



Introduction

Human Alpha-1 Antichymotrypsin (SERPINA3) rabbit polyclonal antibody, suitable for WB, IHC-P, ELISA applications. A secreted serpin implicated in inflammation and Alzheimer's amyloid plaques.

[Learn More](#)

Parameter	Value
Product Name	Alpha-1 Antichymotrypsin (SERPINA3) Rabbit pAb
Remarks/Alias	ACT, AACT, GIG24, GIG25, Alpha-1 Antichymotrypsin (SERPINA3)
Species	Human
Gene ID	12
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 24-423 of human Alpha-1 Antichymotrypsin (SERPINA3) (NP_001076.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P01011
Protein Weight	48kDa
Shipping	Ice bag