

Alpha-1 Antitrypsin (Serpina1) Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P01009

Item Number: CM0018466



Introduction

Rabbit polyclonal antibody against human Alpha-1 Antitrypsin (SERPINA1), suitable for WB, IF/ICC, ELISA. Validated in human, mouse, and rat. Target: serine protease inhibitor, involved in lung protection and hemostasis. UniProt: P01009.

[Learn More](#)

Parameter	Value
Product Name	Alpha-1 Antitrypsin (SERPINA1) Rabbit pAb
Remarks/Alias	PI; A1A; AAT; PI1; A1AT; nNIF; PRO2275; alpha1AT; Alpha-1 Antitrypsin (SERPINA1)
Species	Human
Gene ID (Human)	5265
Gene ID	5265
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-315 of human Alpha-1 Antitrypsin (NP_000286.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01009
Protein Weight	47kDa
Shipping	Ice bag