

Alpha-Fetoprotein (Afp) Rabbit Pab - P02771

Item Number: CM0016278



Introduction

High-quality rabbit polyclonal antibody against Alpha-Fetoprotein (AFP). Validated for WB, IF/ICC, and ELISA with cross-reactivity to human and mouse AFP. Ideal for liver cancer biomarker research and prenatal diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	Alpha-Fetoprotein (AFP) Rabbit pAb
Remarks/Alias	AFPD; FETA; HPAFP; Alpha-Fetoprotein (AFP)
Species	Human
GeneID (Human)	174
GeneID	174
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 450-550 of Alpha-Fetoprotein (AFP) (NP_001125.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P02771
Protein Weight	69kDa
Shipping	Ice bag