

Alpha2-Adaptin/AP2A2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O94973

Item Number: CM0025281

Introduction

Rabbit polyclonal antibody against human Alpha2-adaptin/AP2A2 (AP-2 complex subunit alpha-2). Validated for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Protein weight 104 kDa. Ideal for endocytosis and protein trafficking studies.

[Learn More](#)

Parameter	Specification
Product Name	Alpha2-adaptin/AP2A2 Rabbit pAb
Remarks/Alias	HIP9; HYPJ; ADTAB; HIP-9; CLAPA2; Alpha2-adaptin/AP2A2
Species	Human
GeneID (Human)	161
GeneID	161
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 590-939 of human Alpha2-adaptin/AP2A2 (NP_036437.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O94973
Protein Weight	104kDa
Shipping	Ice bag