

Als2 Rabbit Pab Polyclonal Antibody For Als2 (Alsin) With Applications In Wb, If/Icc, If-P, Elisa - Q96Q42

Item Number: CM0029479

Introduction

Rabbit polyclonal antibody against human ALS2 (alsin) protein, validated for WB, IF/ICC, IF-P, ELISA. Targets an amyotrophic lateral sclerosis-linked protein with cross-reactivity to mouse and rat.

[Learn More](#)

Parameter	Value
Product Name	ALS2 Rabbit pAb
Remarks/Alias	ALSJ; PLSJ; IAHSP; ALS2CR6; ALS2
Species	Human
GeneID (Human)	57679
GeneID	57679
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 260-429 of human ALS2 (NP_065970.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96Q42
Protein Weight	184kDa
Shipping	Ice bag