

Ampd1 Rabbit Pab For Wb, Elisa - P23109

Item Number: CM0030819

Introduction

AMPD1 Rabbit pAb is a polyclonal antibody raised against a synthetic peptide from human AMP deaminase 1. Validated for WB and ELISA, it detects human, mouse, and rat AMPD1, a key enzyme in energy metabolism.

[Learn More](#)

Parameter	Value
Product Name	AMPD1 Rabbit pAb
Remarks/Alias	MAD, MADA, MMDD, AMPD1
Species	Human
GeneID (human)	270
GeneID	270
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human AMPD1 (P23109)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23109
Protein Weight	86kDa
Shipping	Ice bag