

Ampd3 Rabbit Pab For Wb - Q01432

Item Number: CM0025018

Introduction

AMPD3 Rabbit Polyclonal Antibody validated for Western Blot and ELISA. Detects human, mouse, and rat AMP deaminase 3, a key enzyme in energy metabolism. Ideal for nucleotide metabolism research.

[Learn More](#)

Parameter	Value
Product Name	AMPD3 Rabbit pAb
Remarks/Alias	AMPD3
Species	Human
Gene ID (Human)	272
Gene ID	272
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human AMPD3 (NP_000471.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q01432
Protein Weight	89kDa
Shipping	Ice bag