


Anapc10 Rabbit Pab - Q9Um13

Item Number: CM0030647



Introduction

ANAPC10 Rabbit pAb is a polyclonal antibody reactive against human, mouse, and rat ANAPC10 (21kDa), suitable for WB and ELISA. It targets the APC/C subunit 10, a key regulator of cell cycle progression and ubiquitination.

[Learn More](#)

Parameter	Value
Product Name	ANAPC10 Rabbit pAb
Remarks/Alias	DOC1; APC10; ANAPC10
Species	Human
GeneID (Human)	10393
GeneID	10393
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-185 of human ANAPC10 (NP_055700.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UM13
Protein Weight	21kDa
Shipping	Ice bag