

Angel2 Rabbit Polyclonal Antibody For Wb, Elisa - Q5Vte6

Item Number: CM0016265



Introduction

High-quality rabbit polyclonal antibody targeting human ANGEL2 (RNA 2',3'-cyclic phosphatase) validated for Western blot and ELISA. Ideal for mitochondrial RNA processing and tRNA recycling studies.

[Learn More](#)

Parameter	Specification
Product Name	ANGEL2 Rabbit pAb
Remarks/Alias	Ccr4d; KIAA0759L; ANGEL2
Species	Human
Gene ID (Human)	90806
Gene ID	90806
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-170 of human ANGEL2 (NP_653168.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q5VTE6
Protein Weight	62kDa
Shipping	Ice bag