

Annexin A4/Annexin Iv Rabbit Polyclonal Antibody For Wb, Elisa - P09525

Item Number: CM0030268



Introduction

High-quality rabbit polyclonal antibody targeting Annexin A4, validated for WB and ELISA in human, mouse, and rat samples. Useful for studying membrane fusion and exocytosis.

[Learn More](#)

Parameter	Value
Product Name	Annexin A4/Annexin IV Rabbit pAb
Remarks/Alias	ANX4; P32.5; PIG28; PP4-X; ZAP36; PAP-II; HEL-S-274; Annexin A4/Annexin IV
Species	Human
GeneID (Human)	307
GeneID	307
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 112-321 of human Annexin A4/Annexin IV (NP_001144.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09525
Protein Weight	36kDa
Shipping	Ice bag