


Annexin A6/Annexin VI Rabbit Polyclonal Antibody For Wb, Elisa - P08133

Item Number: CM0025202



Introduction

Rabbit polyclonal antibody against human Annexin A6/Annexin VI, validated for WB and ELISA. Cross-reacts with mouse and rat. Suitable for investigating calcium signaling and membrane-related processes.

[Learn More](#)

Parameter	Value
Product Name	Annexin A6/Annexin VI Rabbit pAb
Remarks/Alias	p68; p70; ANX6; CBP68; CPB-II; Annexin A6/Annexin VI
Species	Human
Gene ID	309
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human Annexin A6/Annexin VI (NP_001146.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08133
Protein Weight	76kDa
Shipping	Ice bag