

Anprb/Nrp2 Rabbit Pab - P20594

Item Number: CM0021008

Introduction

Rabbit polyclonal antibody targeting ANPRB/NRP2, validated for WB and ELISA in human samples. This receptor binds CNP and regulates skeletal growth, with a molecular weight of 117 kDa.

[Learn More](#)

Parameter	Value
Product Name	ANPRB/NRP2 Rabbit pAb
Remarks/Alias	GCB; AMD1; AMDM; ANPb; ECDM; GC-B; NPRB; SNSK; ANPRB; GUC2B; NPRBi; GUCY2B; ANPRB/NRP2
Species	Human
Gene ID (Human)	4882
Gene ID	4882
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 130-320 of human ANPRB/NRP2 (NP_003986.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P20594
Protein Weight	117kDa
Shipping	Ice bag