

Ant2/Slc25A5 Rabbit Pab For Wb, If/Icc, Ihc-P, Elisa - P05141

Item Number: CM0017385

Introduction

Rabbit polyclonal antibody against ANT2/SLC25A5, validated for WB, IF/ICC, IHC-P, and ELISA. Recognizes human, mouse, rat, and monkey ANT2. Ideal for mitochondrial apoptosis, metabolism, and mitophagy research.

[Learn More](#)

Parameter	Value
Product Name	ANT2/SLC25A5 Rabbit pAb
Remarks/Alias	T2; T3; 2F1; AAC2; ANT2; SLC25A5
Species	Human
GeneID (Human)	292
GeneID	292
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat, Monkey
SWISS	P05141
Protein Weight	33 kDa
Shipping	Ice bag