

Anti-14-3-3 Alpha/Beta Monoclonal Antibody For Wb, Elisa - P31946

Item Number: CM0010492



Introduction

Rabbit monoclonal antibody targeting human 14-3-3 alpha/beta (YWHAB). Validated in WB and ELISA; cross-reacts with mouse and rat. Ideal for apoptosis, osteogenesis, and MAP kinase signaling studies.

[Learn More](#)

Parameter	Value
Product Name	14-3-3 alpha/beta Rabbit mAb
Remarks/Alias	HS1; GW128; YWHAA; KCIP-1; HEL-S-1; 14-3-3 alpha/beta
Species	Human
GeneID	7529
Immunogen	Synthetic Peptide - A synthetic peptide corresponding to a sequence within amino acids 147-246 of human 14-3-3 alpha/beta (P31946).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P31946
Protein Weight	28kDa
Shipping	Ice bag