


Anti-14-3-3 Alpha/Beta Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - P31946

Item Number: CM0016163



Introduction

Anti-14-3-3 alpha/beta Rabbit Polyclonal Antibody validated for WB, IHC-P, IF/ICC, IP, ELISA. Targets human 14-3-3 beta/alpha (P31946), with cross-reactivity to mouse and rat. Ideal for signaling pathway research.

[Learn More](#)

Parameter	Value
Product Name	14-3-3 alpha/beta Rabbit pAb
Remarks/Alias	HS1; GW128; YWHAA; KCIP-1; HEL-S-1; 14-3-3 alpha/beta
Species	Human
GeneID (Human)	7529
GeneID	7529
Immunogen	Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 150-246 of human 14-3-3 alpha/beta (NP_003395.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P31946
Protein Weight	28kDa
Shipping	Ice bag