

Anti-14-3-3-Beta/Zeta Rabbit Monoclonal Antibody For Wb, Elisa - P31946 / P63104

Item Number: CM0004100



Introduction

High-quality rabbit monoclonal antibody against human 14-3-3 beta/zeta (YWHAB/YWHAZ). Validated for WB and ELISA, cross-reacts with mouse and rat. Ideal for signaling pathway research.

[Learn More](#)

Parameter	Value
Product Name	14-3-3-beta/zeta Rabbit mAb
Remarks/Alias	GW128; HEL-S-1; HS1; KCIP-1; YWHAA; 14-3-3-beta/zeta
Species	Human
GeneID (Human)	7529/7534
Gene ID	7529/7534
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human 14-3-3 beta/zeta (P31946).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P31946/P63104
Protein Weight	28kDa
Shipping	Ice bag