

Anti-14-3-3 Beta/Zeta Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P31946 / P63104

Item Number: CM0009401



Introduction

Rabbit monoclonal antibody against human 14-3-3 beta/zeta (P31946/P63104), validated for WB, IF/ICC, ELISA. Recognizes 28 kDa protein; cross-reacts with mouse and rat. Suitable for signaling pathway research.

[Learn More](#)

Parameter	Value
Product Name	14-3-3 beta/zeta Rabbit mAb
Remarks/Alias	GW128; HEL-S-1; HS1; KCIP-1; YWHAA; 14-3-3; 14-3-3 beta/zeta
Species	Human
GeneID (human)	7529/7534
GeneID	7529/7534
Immunogen	Synthetic Peptide (A synthetic peptide corresponding to a sequence within amino acids 1-100 of human 14-3-3 beta/zeta (P31946))
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P31946/P63104
Protein Weight	28kDa
Shipping	Ice bag