


Anti-14-3-3 Epsilon Monoclonal Antibody For Wb, Elisa - P62258

Item Number: CM0008120



Introduction

Rabbit monoclonal antibody against human 14-3-3 epsilon (P62258), validated for WB and ELISA. Suitable for detecting endogenous 14-3-3 epsilon in human, mouse, and rat samples.

[Learn More](#)

Parameter	Value
Product Name	14-3-3 epsilon Rabbit mAb
Remarks/Alias	MDS; HEL2; MDCR; KCIP-1; 14-3-3E; 14-3-3 epsilon
Species	Human
GeneID (human)	7531
GeneID	7531
Immunogen	Synthetic Peptide - A synthetic peptide corresponding to a sequence within amino acids 1-100 of human 14-3-3 epsilon (P62258)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62258
Protein Weight	28kDa
Shipping	Ice bag