


Anti-14-3-3 Epsilon Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P62258

Item Number: CM0016427



Introduction

Rabbit polyclonal antibody against human 14-3-3 epsilon (29kDa), validated in WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Recombinant immunogen (aa 73-163). Supports research into cell signaling, antiviral pathways, and protein interactions.

[Learn More](#)

Parameter	Value
Product Name	14-3-3 epsilon Rabbit pAb
Remarks/Alias	MDS; HEL2; MDCR; KCIP-1; 14-3-3E; 14-3-3 epsilon
Species	Human
GeneID (Human)	7531
GeneID	7531
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 73-163 of human 14-3-3 epsilon (NP_006752.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62258
Protein Weight	29kDa
Shipping	Ice bag