


Anti-14-3-3 Gamma Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P61981

Item Number: CM0011085



Introduction

Rabbit monoclonal antibody against human 14-3-3 gamma (P61981). Validated for WB, IHC-P, ELISA with cross-reactivity to human, mouse, rat. Target molecular weight 28kDa. Ideal for signaling pathway research.

[Learn More](#)

Parameter	Value
Product Name	14-3-3 gamma Rabbit mAb
Remarks/Alias	DEE56; EIEE56; PPP1R170; 14-3-3GAMMA; 14-3-3 gamma
Species	Human
GeneID (Human)	7532
GeneID	7532
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 100-200 of human 14-3-3 gamma (P61981).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P61981
Protein Weight	28kDa
Shipping	Ice bag