


## Anti-14-3-3 Sigma Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P31947

Item Number: CM0005375



### Introduction

High-quality anti-14-3-3 sigma rabbit monoclonal antibody validated for Western blot, IHC-P, and ELISA. High specificity for human and mouse 14-3-3 sigma (P31947), an adapter protein regulating cell cycle and signaling.

[Learn More](#)

| Parameter        | Value                                                                                                                      |
|------------------|----------------------------------------------------------------------------------------------------------------------------|
| Product Name     | 14-3-3 sigma Rabbit mAb                                                                                                    |
| Remarks/Alias    | YWHAS; 14-3-3 sigma                                                                                                        |
| Species          | Human                                                                                                                      |
| GeneID (Human)   | 2810                                                                                                                       |
| GeneID           | 2810                                                                                                                       |
| Immunogen        | Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human 14-3-3 sigma (P31947). |
| Source           | Rabbit                                                                                                                     |
| Category         | Monoclonal Antibodies                                                                                                      |
| Application      | WB, IHC-P, ELISA                                                                                                           |
| Cross Reactivity | Human, Mouse                                                                                                               |
| SWISS            | P31947                                                                                                                     |
| Protein Weight   | 28kDa                                                                                                                      |
| Shipping         | Ice bag                                                                                                                    |