


Anti-14-3-3 Sigma Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P31947

Item Number: CM0018292



Introduction

Rabbit polyclonal antibody targeting human 14-3-3 sigma (stratifin) for WB, IHC-P, and ELISA. Cross-reacts with mouse and rat. Ideal for studying cell signaling, cell cycle, and epithelial biology.

[Learn More](#)

Parameter	Value
Product Name	14-3-3 sigma Rabbit pAb
Remarks/Alias	YWHA5; 14-3-3 sigma
Species	Human
Gene ID	2810
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human 14-3-3 sigma (NP_006133.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P31947
Protein Weight	28kDa
Shipping	Ice bag