


Anti-14-3-3 Theta Polyclonal Antibody For Wb, Elisa - P27348

Item Number: CM0022464



Introduction

Rabbit polyclonal antibody targeting human 14-3-3 theta (YWHAQ) for Western blot and ELISA applications. Cross-reacts with mouse and rat. 28 kDa signaling adapter protein. P27348.

[Learn More](#)

Parameter	Value
Product Name	14-3-3 Theta Rabbit pAb
Remarks/Alias	1C5; HS1; 14-3-3; 14-3-3 Theta
Species	Human
GeneID	10971
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 49-149 of human 14-3-3 Theta (NP_006817.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P27348
Protein Weight	28kDa
Shipping	Ice bag