


Anti-14-3-3 Theta Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P27348

Item Number: CM0011065



Introduction

High-quality rabbit monoclonal antibody targeting 14-3-3 theta (YWHAQ). Validated for Western blot, IF/ICC, and ELISA with cross-reactivity to human, mouse, and rat. Ideal for cell signaling and neurodegenerative disease research.

[Learn More](#)

Parameter	Specification
Product Name	14-3-3 Theta Rabbit mAb
Remarks/Alias	1C5; HS1; 14-3-3; 14-3-3 Theta
Species	Human
Gene ID (Human)	10971
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human 14-3-3 Theta (P27348)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P27348
Protein Weight	28kDa
Shipping	Ice bag