


Anti-14-3-3 Zeta/Ywhaz Monoclonal Antibody For Wb, If/Icc, Elisa - P63104

Item Number: CM0007426



Introduction

Rabbit monoclonal antibody targeting 14-3-3 zeta/YWHAZ (P63104). Suitable for WB, IF/ICC, ELISA applications, with cross-reactivity to human, mouse, and rat. Ideal for studying adapter protein signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	14-3-3 zeta/YWHAZ Rabbit mAb
Remarks/Alias	HEL4; YWHAD; KCIP-1; HEL-S-3; POPCHAS; HEL-S-93; 14-3-3-zeta; 14-3-3 zeta/YWHAZ
Species	Human
GeneID	7534
Immunogen	Synthetic Peptide - A synthetic peptide corresponding to a sequence within amino acids 145-245 of human YWHAZ (NP_002267.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P63104
Protein Weight	19kDa/28kDa
Shipping	Ice bag