


Anti-15-Pgdh/Hpgd Polyclonal Antibody For Wb, Elisa - P15428

Item Number: CM0031038



Introduction

Rabbit polyclonal antibody against human 15-hydroxyprostaglandin dehydrogenase, validated in WB and ELISA, cross-reacts with mouse and rat. Suitable for research on prostaglandin metabolism and inflammation.

[Learn More](#)

Parameter	Value
Product Name	15-PGDH/HPGD Rabbit pAb
Remarks/Alias	PGDH; PGDH1; PHOAR1; 15-PGDH; SDR36C1; 15-PGDH/HPGD
Species	Human
GeneID	3248
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-266 of human 15-PGDH/HPGD (NP_000851.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P15428
Protein Weight	29kDa
Shipping	Ice bag