

Anti-4-1Bb/Cd137 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q07011

Item Number: CM0023006



Introduction

This rabbit polyclonal antibody targets human 4-1BB/CD137 (TNFRSF9), a T-cell costimulatory receptor, for applications including WB, IHC-P, and ELISA. Cross-reacts with human, mouse, and rat. Enhances CD8+ T-cell survival and anti-tumor immunity.

[Learn More](#)

Parameter	Value
Product Name	4-1BB/CD137 Rabbit pAb
Remarks/Alias	ILA; 4-1BB; CD137; CDw137; IMD109; 4-1BB/CD137
Species	Human
GeneID (Human)	3604
GeneID	3604
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 100-200 of human 4-1BB/CD137 (NP_001552.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q07011
Protein Weight	28kDa
Shipping	Ice bag