

Anti-4-1Bbl/Cd137L Rabbit Monoclonal Antibody For Wb, Elisa - P41273

Item Number: CM0010635



Introduction

Rabbit monoclonal antibody specific for human 4-1BBL/CD137L (TNFSF9). Validated for WB and ELISA, cross-reacts with mouse and rat. Recognizes a 27kDa protein. Useful for T-cell biology and immuno-oncology studies.

[Learn More](#)

Parameter	Value
Product Name	4-1BBL/CD137L Rabbit mAb
Remarks/Alias	CD137L; TNLG5A; 4-1BB-L; 4-1BBL/CD137L
Species	Human
GeneID (Human)	8744
GeneID	8744
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 155-254 of human 4-1BBL/CD137L (P41273).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P41273
Protein Weight	27kDa
Shipping	Ice bag