


Anti-53BP1 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q12888

Item Number: CM0025982



Introduction

Rabbit polyclonal antibody against human 53BP1 (TP53BP1), reactive with human, mouse, rat. Suitable for WB, IHC-P, ELISA. Key role in DNA damage response and repair, telomere maintenance, and class-switch recombination.

[Learn More](#)

Parameter	Value
Product Name	53BP1 Rabbit pAb
Remarks/Alias	p202; 53BP1; TDRD30; p53BP1
Species	Human
Gene ID (Human)	7158
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 895-1031 of human TP53BP1 (NP_005648.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q12888
Protein Weight	214kDa
Shipping	Ice bag