

Anti-53BP1 Rabbit Monoclonal Antibody For IHC-P/IF/ICC/Elisa - Q12888

Item Number: CM0004439



Introduction

Rabbit monoclonal antibody against human 53BP1, validated for IHC-P, IF/ICC and ELISA. Detects a 214kDa protein; cross-reacts with mouse. Suitable for DNA damage and repair research.

[Learn More](#)

Parameter	Value
Product Name	53BP1 Rabbit mAb
Remarks/Alias	p202; 53BP1; TDRD30; p53BP1
Species	Human
Gene ID	7158 (Human)
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human 53BP1 (Q12888)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q12888
Protein Weight	214kDa
Shipping	Ice bag