

Anti-5Ht7 Receptor Rabbit Monoclonal Antibody For Wb, If-P, Elisa - P34969

Item Number: CM0002396



Introduction

High-specificity rabbit monoclonal antibody targeting human 5HT7 receptor (5-HT7, serotonin receptor 7). Validated for WB, IF-P, ELISA; cross-reacts with mouse and rat. Ideal for investigating serotonin-mediated GPCR signaling.

[Learn More](#)

Parameter	Value
Product Name	5HT7 Receptor Rabbit mAb
Remarks/Alias	5-HT7; 5HT7 Receptor
Species	Human
GeneID (Human)	3363
GeneID	3363
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human 5HT7 Receptor (P34969).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P34969
Protein Weight	54kDa
Shipping	Ice bag