

Anti-5T4 (Trophoblast Glycoprotein) Polyclonal Antibody For Wb, Elisa - Q13641

Item Number: CM0018039



Introduction

Rabbit polyclonal antibody targeting human 5T4 (trophoblast glycoprotein). Validated for WB and ELISA applications with cross-reactivity to mouse and rat. Recognizes the 46 kDa Wnt signaling inhibitor for research use.

[Learn More](#)

Parameter	Value
Product Name	5T4 Rabbit pAb
Remarks/Alias	5T4; M6P1; 5T4AG; WAIF1
Species	Human
GeneID	7162
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 32-350 of human 5T4 (NP_006661.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q13641
Protein Weight	46kDa
Shipping	Ice bag