


Anti-67Kda Laminin Receptor Polyclonal Antibody For Wb, If/Icc, Elisa - P08865

Item Number: CM0028435



Introduction

Anti-67kDa Laminin Receptor Polyclonal Antibody validated in WB, IF/ICC, ELISA with human, mouse, rat reactivity. Explore 67LR's roles in ribosomal assembly, cell adhesion, and viral infection.

[Learn More](#)

Parameter	Value
Product Name	67kDa Laminin Receptor Rabbit pAb
Remarks/Alias	SA; LBP; LRP; p40; uS2; 67LR; ICAS; lamR; 37LRP; LAMBR; LAMR1; LRP/LR; LBP/p40; NEM/1CHD4; 67kDa Laminin Receptor
Species	Human
GeneID (Human)	3921
GeneID	3921
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-295 of human 67kDa Laminin Receptor (NP_002286.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08865
Protein Weight	33kDa
Shipping	Ice bag