


Anti-67Kda Laminin Receptor (Rpsa) Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P08865

Item Number: CM0012473



Introduction

Rabbit monoclonal antibody targeting 67kDa laminin receptor (RPSA, 37LRP) for WB, IF/ICC, ELISA. Recognizes human, mouse, rat. Immunogen: synthetic peptide within aa 196-295. Ideal for cancer and infectious disease research.

[Learn More](#)

Parameter	Value
Product Name	67kDa Laminin Receptor Rabbit mAb
Remarks/Alias	SA; LBP; LRP; p40; uS2; 67LR; ICAS; lamR; 37LRP; LAMBR; LAMR1; LRP/LR; LBP/p40; NEM/1CHD4; 67kDa Laminin Receptor
Species	Human
Gene ID	3921
Immunogen	Synthetic Peptide - A synthetic peptide corresponding to a sequence within amino acids 196-295 of human Laminin Receptor (P08865).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08865
Protein Weight	33kDa
Shipping	Ice bag