

Anti-A-Actin-1 (A-Sma) Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P68133

Item Number: CM0004393



Introduction

Rabbit monoclonal antibody against human α -Actin-1 (α -SMA), validated for WB, IF-P, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for studying muscle biology and cytoskeleton research.

[Learn More](#)

Parameter	Value
Product Name	α -Actin-1 (α -SMA) Rabbit mAb
Remarks/Alias	ACTA; ACTA1; CFTD; MPFD; NEM1; NEM2; NEM3; SHPM; CFTD1; CFTDM; CMYP2A; CMYP2B; CMYP2C; α -Actin-1 (α -SMA)
Species	Human
GeneID (Human)	58
GeneID	58
Immunogen	Synthetic peptide corresponding to amino acids 278-377 of human α -Actin-1 (ACTA1) (P68133)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P68133
Protein Weight	42kDa
Shipping	Ice bag