

Anti-A-Actin-1 (Acta1) Rabbit Monoclonal Antibody For Wb, If/Icc, If-P, Elisa - P68133

Item Number: CM0007180



Introduction

High-quality rabbit monoclonal antibody targeting human α -Actin-1 (ACTA1) for Western blotting (WB), immunofluorescence (IF/ICC, IF-P), and ELISA. Cross-reacts with mouse and rat. Essential tool for studying skeletal muscle cytoskeleton and cell motility.

[Learn More](#)

Parameter	Value
Product Name	α -Actin-1 (ACTA1) Rabbit mAb
Remarks/Alias	ACTA; ASMA; CFTD; MPFD; NEM1; NEM2; NEM3; SHPM; CFTD1; CFTDM; CMYP2A; CMYP2B; CMYP2C; α -Actin-1 (ACTA1)
Species	Human
Gene ID (Human)	58
Gene ID	58
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human α -Actin-1 (ACTA1) (P68133)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P68133
Protein Weight	42kDa
Shipping	Ice bag