


Anti-A-Actinin-4 Polyclonal Antibody For Wb, Elisa - O43707

Item Number: CM0023323



Introduction

Rabbit polyclonal antibody against human α -Actinin-4, validated for WB and ELISA. Cross-reacts with human, mouse, and rat. Suitable for research on F-actin cross-linking, vesicular trafficking, and tight junction biology.

[Learn More](#)

Parameter	Value
Product Name	α -Actinin-4 Rabbit pAb
Remarks/Alias	FSGS; FSGS1; ACTININ-4; α -Actinin-4
Species	Human
GeneID (Human)	81
GeneID	81
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 570-650 of human α -Actinin-4 (NP_004915.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43707
Protein Weight	105kDa
Shipping	Ice bag