

Anti-A-Raf Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P10398

Item Number: CM0005584



Introduction

Rabbit monoclonal antibody against human A-RAF (P10398). Validated for Western blot (WB), IHC-P, and ELISA. Suitable for studying mitogenic signaling, myogenic differentiation, and cancer pathways.

[Learn More](#)

Parameter	Value
Product Name	A-RAF Rabbit mAb
Remarks/Alias	PKS2; A-RAF; ARAF1; RAFA1
Species	Human
GeneID (Human)	369
GeneID	369
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 166-290 of human A-RAF (P10398)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P10398
Protein Weight	68kDa
Shipping	Ice bag