

Anti-A-Smooth Muscle Actin (Acta2) Rabbit Monoclonal Antibody For Wb, Elisa - P62736

Item Number: CM0002732

Introduction

Rabbit monoclonal antibody against human α -Smooth Muscle Actin (ACTA2), validated for WB and ELISA. Cross-reacts with mouse and rat. Suitable for cytoskeleton and smooth muscle research.

[Learn More](#)

Parameter	Value
Product Name	α -Smooth Muscle Actin (ACTA2) Rabbit mAb
Remarks/Alias	ACTSA; α -Smooth Muscle Actin (ACTA2)
Species	Human
GeneID (Human)	59
GeneID	59
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human α -Smooth Muscle Actin (ACTA2) (P62736)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62736
Protein Weight	42kDa
Shipping	Ice bag